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
Understanding the function mode of the two sides of Eye/U

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“**Paste it on where it hurts**” is the basic principle of the use of the Eye/ U, it could be adopted the methods of “Paste two Eye/U on two against body sides” and “Paste two Eye/U with energy-in and one with energy-out”.

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
“Paste two Eye/U on two against body sides” means that paste two Eye/U with both energy-in or energy-out on the against body sides around the pain area. Function A

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“Paste two Eye/U with energy-in and one with energy-out” means that paste two Eye/U with energy-in on the near-end around far from the painful area 1 centimeter of the limbs, paste one Eye/U with energy-out on the far-end around far from the painful area 1 centimeter of the limbs. Function B

Conditioning to knee and leg pain with Function A & B

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How to find out the main trunk line of the human-channel and the acupoints precisely?

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The main trunk lines of the human-channel inherently distribute along the edge of the limbs bone. Their distribution in head, neck and torso is complex. The acupoints distribute all over the body which commonly present a tissue morphology with a thin center and thick periphery.

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Recommended conditioning schemes

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
Conditioning to large intestine related diseases as an Assistant

Conditioning to large intestine related diseases as an Assistant

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
Diseases related to the large intestine include Hemorrhoids, Crohn's disease, Constipation, Enteritis, Intestinal obstruction and colon cancer. Because there is a close relationship between the large intestine and the lungs in terms of Human-channel(meridian) connection, when treating lung diseases, it is generally advisable to also take into account the regulation of the large intestine. The conditioning method shown in the above picture is a general one. To achieve the best conditioning effect, other auxiliary conditioning plans are still needed. For technical support, please contact us.

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
Conditioning to heart discomfort after vaccination

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Conditioning to the neck or head pain


Conditioning to the neck or head pain

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Enhanced Conditioning to facial paralysis as an Assistant


Enhanced Conditioning to facial paralysis as an Assistant

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Conditioning for pain and itch relief

Conditioning for pain and itch relief

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
The first step is to paste the Eye/U to the position that needs to be adjusted with a transparent adhesive patch positioning. The second step is to cover the transparent adhesive patch with an elastic fixing patch for reinforcement.

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Request technical support is the best way to have a good experience. The functions of the human body always show a pattern of being sometimes strong and then weak, and the symptoms of diseases exhibit characteristics of being mild, severe and relieved. These patterns are

caused by individual organs being on duty of president during the operation of the Organ Rhythm. Seizing the Handover Period, the Preparing Period and the President Period of the relevant organs to carry out targeted conditioning is a highly efficient method to improve the operation efficiency of the Organ Rhythm and self-healing ability. Under normal circumstances, a period of professional technical support may be helpful to overcome difficult diseases, or it may be helpful to achieve beauty and body shaping by “stimulating intrinsic motivation of the human body”. Please pay attention to our website, we will publish recommended conditioning schemes from time to time, or contact us to learn about the conditioning value-added project.

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**The Eye/U and the Organ Rhythm Conditioner
provides a strong health support to your family, Effect Talks!**

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Beware of adhesive patch allergy! If you are allergic to the adhesive patch, you can use gel to stick the Eye/U onto your skin.

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Beijing Wanhua Technology Co., Ltd