



Organ Rhythm Conditioner



Instruction
manual

United States of America

United States Patent and Trademark Office



Reg. No. 7,252,348

Registered Dec. 26, 2023

Int. Cl.: 10

Trademark

Principal Register

Beijing Zhihuishan Technology Co., Ltd. (CHINA limited company (ltd.))

Rm 707, 7th Floor, Building 689, Dist 1,

No. 5 South Street, Haidian Dist,

Beijing, CHINA 100089

CLASS 10: Acupuncture instruments; Acupuncture needles; Apparatus for acupuncture therapy; Apparatus for magnet therapy; Body rehabilitation apparatus for medical purposes; Electromedical rehabilitative and pain management products for clinical and home use, namely, electrical nerve and muscle stimulators, ultrasonic stimulators, magnet therapy stimulators and laser therapy stimulators; Low frequency electric therapy apparatus; Low frequency electromagnetic therapy apparatus; Medical and therapeutic device and apparatus, namely, a force and motion sensing apparatus and assisted exercise machine for the rehabilitation of hands and wrists, feet and ankles, knees affected by neuromuscular diseases, disorders, or injuries for home or clinical use; Medical device for accessing subcutaneous structures of the human body for the purpose of administering medical treatment; Medical devices for nonsurgical cosmetic treatments; Moxibustion apparatus; Non-medicated compresses; Portable medical devices with sensors to monitor the physical movements of a patient wearing or carrying the device; Silicone gel sheeting for the treatment of scars; Therapeutic weighted products, namely, weighted blankets, lap pads and vests designed for clinical and home use to provide proprioceptive input to the body

FIRST USE 4-26-2021; IN COMMERCE 4-26-2021

The mark consists of stylized wording "WAN HUA". On the left is part of the meridian structure pattern in the human small intestine, on the right are two irregular Chinese characters, and the word "Wanhua" under the first Chinese character.

The English translation of "WAN HUA" in the mark is "ten thousand" and "digest".

The non-Latin characters in the mark transliterate to "wan" and "hua" and this means "ten thousand" and "digestion" in English.

Katherine Kelly Vidal

Director of the United States
Patent and Trademark Office



CATALOGUE

Introduction	01
Cautions	02
Product diagram	03
1. Host interface of the Organ Rhythm Conditioner	
2. The description of the wrist strap and ankle strap	
Instructions of operation and display	05
How to wear the wrist strap and ankle strap	07
1. The correspondence of the gold-fingers and the human channels	
2.The specific method of the position adjustment of the gold-finger	
A functional description of the conditioning program	11
Specification parameter	12
Declaration	12
After service	12
Frontier of life science	13
The answers to the problems of health	14

INTRODUCTION

The Organ Rhythm Conditioner(ORC) is a health conditioning equipment developed according to the latest research results of the "Rotating Presidency of Twelve-organ", which reveals the physiological rhythm and human channel operation characteristics. This health conditioning method of the ORC is based on the basic principle above, the treatment does not invade into human tissues, without any stimulation. It is a human-friendly wearable instrument for all-weather, and it can be used in many scenarios.

The principle of health regulating is that when the wrist straps and ankle straps are fixed to the wrist and ankle, the gold-fingers are aligned with the twelve channels on the skin surface, the ORC can lead the directional changing magnetic field into the target channels, and correct the flowing direction of JingQi (the gaseous substance flowing in the human channel) which needs to be cleared or promoted back to normal status. It can expand the channels to promote JingQi diffusion well in local tissues and cells through promoting the transverse flow. JingQi in both sides of their channels are in reverse flow state when the twelve organs are on duty of the President period. The synergism of metabolism in various parts of human body can be effectively regulated if the reverse flow of JingQi could be promoted smoothly during the President period because there becomes a temporary circulation around the president organ. That is why the optimal human function of each organ with its own characteristics can be enhanced.

In the process of human body from birth to aging, it is not possible to completely shield the invasion of wind, cold, hot, dampness, dry and fire of the environment, and it is difficult to ensure that the ingested grains and cereals are fully digested and perfectly utilized by the body, then the resulting deviation in the rhythm of organ and the operation of the channels will inevitably be manifested in the way of disease. Common hypertension, heart disease, diabetes, sleep disorders, and all kinds of difficult and complicated diseases are the reflection of the disturbance of the organ rhythm and blockage of the channels. Therefore, the daily use of the ORC for organ rhythm and human channels conditioning can effectively ensure the normal operation of human growth, digestion, collection and immune self-healing function, so as to achieve the self-repair of the human's sub-health and diseased states.

It can be said that the body has been relaxed and feeling well with having no time to sensing.

Although the Organ Rhythm Conditioner has a very excellent role in health maintenance, but it is not a substitute for medical treatment, please remember!


CAUTIONS

Please be sure to use the charging head with output voltage no higher than DC5V, and do not use the fast charger of no higher than 2A. Otherwise it may cause the machine overheat to burn.

The company reserves the rights to improve the product, specifications and design of the product, any subject changes will not notice. The pictures and text descriptions in this manual are for reference only, subject to the application requirements is as the actual use of the product.

PRODUCT DIAGRAM

1. Host interface of the Organ Rhythm Conditioner



1. Display screen

2. Power button

3. Function button

4. Connection Cables


5. Memory card slot

6. Rotatable back clip


7. Product parameter

8. Charging port


9. Charging cable



2. The description of the wrist strap and ankle strap



Outside of the wristband and the ankle



Inside of the wristband and the ankle



- 1. TYPE C socket
- 2. Iron hoop
- 3. HOOK&LOOP

- 4. Cable tie
- 5. Gold-finger
- 6. Retaining clip


INSTRUCTIONS OF OPERATION AND DISPLAY

Power on




Press and hold the "Power" button until the display lights up and releases, about 3 seconds. It displays “Healthy me Ready” and battery icon.


Program conversion



① Click the "Function" button, the program number is switched in ascending order. For example, click the "Function" button to display program 1, and there is a dynamic display of the usage time in program 1; Click the "Function" button again to display program 2. After reaching the highest program number in the equipment, click again to jump back to the interface of “Healthy me Ready” and battery icon.



② When an even number of program is displayed, there is a dynamic display of the usage time in the program and a motor rotation sign.



③ The program number you switch by clicking the "Function" button is related to the number of programs you have purchased, and the program number is not always in continuous increments. Double tap the "Function" button consecutively to return to program 1.

Close the program



Double tap the "Power" button consecutively to end the running program, the display will display the interface of "Healthy me Ready" and battery icon.


Turning off the screen

Keep pressing the "Power" button until the display goes off. At this time, the program continues to execute, but the motor rotation stops, and the LED inside the equipment lights up.

Wake up the resting screen

keep pressing the "Power" button until the display continues to light up and displays the program number.

Power Off




Press the "Function" and "Power" buttons at the same time until the text of "power off" appears.

HOW TO WEAR THE WRIST STRAP AND ANKLE STRAP



1.The correspondence of the gold-fingers and the human channels

Wrist and ankle straps must be worn with the gold face of the gold-finger aligned to the human channels. The channels are located at the bone edges of the wrist and ankle, as shown below:



In Figure 1- Figure 2, the sequence of the golden-finger from 1-6 aligning to the channels on wrist is as follows:

golden-finger 1- Heart Channel, golden-finger 2- Pericardium Channel, gold-finger 3- Lung Channel, gold-finger 4- Large Intestine Channel, gold-finger 5- Sanjiao Channel, gold-finger 6- Small Intestine Channel.




In Figure 3- Figure 4, the sequence of the golden-finger from 1-6 aligning to the channels on ankle is as follows:

gold-finger 1- Kidney Channel, gold-finger 2- Spleen/Pancreas Channel, gold-finger 3- Liver Channel, gold-finger 4- Stomach Channel, gold-finger 5- Gallbladder Channel, gold-finger 6- Bladder Channel.

2.The specific method of the position adjustment of the gold-finger

2.1 The bottom of the gold-finger has hook that can be used to attach the gold-finger module to the wrist strap and ankle strap, and to adjust the position of the gold-finger on the wrist strap and ankle strap. Be sure to buckle the hook hair at the bottom of the gold-finger with your finger and then pull it, otherwise it is easy to separate the hook from the air bag and cause damage! After separating the gold-finger from the wrist strap and ankle strap, adjust the position of the gold-finger on the wrist strap and ankle strap according to the above figure, so that the gold-fingers can align with each channel after the wrist strap and ankle strap are wrapped around the wrist and ankle nicely, and then the adjustment is completed by putting the extra lead cable into the wrist strap and ankle strap.



2.2 The Organ Rhythm Conditioner is wearable and can be used in various scenarios such as fitness, commuting, home and sleep.




Fitness and commuting use, need to properly adjust the cable, so as not to affect the activity. At this time, you need to tie the cable to the inside of the coat.

2.2.1 First identify the left wrist strap, right wrist strap, left ankle strap, and right ankle strap according to the sign of the TYPE C socket. If you are wearing it by yourself, it is best to insert the round end of the wrist strap into the iron hoop in advance and then put it around the wrist, the TYPE C socket is on the inner wrist and the socket port is facing the body. While pulling the wrist strap tight, check whether the gold face of each gold-finger is properly in contact with the corresponding channels on the skin. A little rotation while pressing on the HOOK&LOOP will make the bond fix stronger.

2.2.2 The first step is to wrap ankle strap around the ankle, then insert the round end of the ankle strap into the iron hoop, the TYPE C socket is located on the inner ankle, and the socket port is facing upward. The following steps are the same as the wrist strap.

2.2.3 Slide the four cable ties to the bottom of the ORC, fasten the free end, and secure the retaining clip into the long hole at the root of the cable ties. Then distinguish the cables of the left hand, right hand, left foot and right foot according to the sign of the cable tie, put the cables throughout the sleeves and trousers respectively, insert the plug into the corresponding socket of the wrist strap and ankle strap, it is properly connected while hearing a click. Finally, clamp the ORC on the waistband or clothes, adjust the position of the cable tie to reserve a certain length of cable for the activities of the arms and legs, and make the retaining clip to the appropriate position of the clothes. At this time, the wearing of the ORC is completed.

2.3 The activities use at home and during sleep are not intense. For the convenience of wearing and taking off, you can place the cable on the outside of the jacket. The steps for adjusting the position of the cable tie and properly fixing it are the same as above. If you wear it by yourself, it is best to connect the plug into the socket on the wrist strap first, because it is difficult to plug in with one hand.



A FUNCTIONAL DESCRIPTION OF THE CONDITIONING PROGRAM

This Organ Rhythm Conditioner has 16 built-in programs. The odd-numbered program is used to improve the coordination of the relevant organs and channels, and to increase energy and the content of JingQi. The even-numbered program is used to improve the coordination of relevant organs and channels, to strengthen the diffusion of JingQi in the whole body and expand the channels.

Prohibited: Without professional guidance, it is forbidden to use even-numbered programs for more than 2 hours continuously!!!

Program 1-2, Improve the coordination of organs and channels.

Program 3-4, Cardio-centered metabolism conditioning.

Program 5-6, Digestion-centered metabolism conditioning.

Program 7-8, Comprehensive conditioning centered on the cardiovascular system.

Program 9-10, Cardiac circuit conditioning.

Program 11-12, Improve the coordination of five-organ and channels.

Program 13-14, Adjust to the unbalanced channels on both sides of the body

Program 15-16, The improper use of these two programs is easy to cause headaches and other uncomfortable symptoms, it must be used under the guidance of professional staff of the manufacturer! These two programs are used for immune pathway dredging, regulation and immune functions strengthening.

SPECIFICATION PARAMETER



DECLARATION

Product packaging accessories and built-in programs are subject to change without prior notice. Please pay attention to the company's website <https://www.wanhua.life> for details of services and offers.

AFTER SERVICE

If the purchased equipment is damaged or the battery is aging during use, please contact our after-sales service agency or dealer, we will solve it as soon as possible.


FRONTIER OF LIFE SCIENCE

Beijing Wanhua Technology Co., Ltd. has deeply researched in human life science for many years. In the comparative study of the Circadian Rhythm of modern medicine and the human channel theory of ZiWuLiuZhu of TCM, the common underlying physiological mechanism to the above two medical systems has been found, which is the rhythm of Rotating Presidency of Twelve-organ(RPTO). The rhythm of RPTO reveals the detailed rotating process of the Handover period, the Preparing period and the President period, as well as the supporting pathway structure of the rhythm operation of RPTO. The findings provide a perfect interpretation for understanding the optimal physiological performance of Circadian Rhythm, and provides a new perspective for understanding health and disease, as well as various life phenomena. In particular, the phenomenon of substances transfer between relevant organ's and the diffusion of the characteristic functional substances to tissues and cells throughout the body during the President period can provide the pathological hints and the healing opportunity.

We specialize in the etiology analysis of complex and rare diseases, and we have developed the association analysis method for multi-factors: Characteristic behavior/ Disease symptoms/ Damaged organs/ Time/ Optimal physiological function/ Characteristic functional substances/ Circulation path/ JingQi flowing direction. Based on this analysis method, many diseases can be three-dimensional images on the paper.

When combined with the use of the product of Eye of the Universe(Eye/U) that activates and opens up the blocked acupoints under our guidance, you could quickly gain improvements in the rhythm of RPTO and the normal status of JingQi flowing. So that make you embrace health through restoring your inherent ability of 'Born, Growth, Digestion, Collection, Reservation and Immunity'.


Welcome to inquire!



Eye of the Universe(Eye/U)

THE ANSWERS TO THE PROBLEMS OF HEALTH

The answers to the problems of your health are right here!



Beijing Wanhua Technology Co., Ltd

Address: Room 1113, 11th Floor, Building 689, District 1,
No.5 Zhongguancun South Street, Haidian District, Beijing

Website: <https://www.wanhua.life>

Email: wanhuatech.sina.com